

Vernier Energy Sensor

(Order Code VES-BTA)



The Vernier Energy Sensor allows students to easily measure current and voltage. Source terminals connect to energy output sources such as model wind turbines or solar panels, and Load terminals connect to loads such as LEDs, water pumps, resistors, or variable loads.

Note: Vernier products are designed for educational use. Our products are not designed nor are they recommended for any industrial, medical, or commercial process such as life support, patient diagnosis, control of a manufacturing process, or industrial testing of any kind.

Compatible Software

See www.vernier.com/manuals/ves-bta for a list of software compatible with the Vernier Energy Sensor.

Quick Start

1. Plug the sensor into the interface (LabQuest 3, LabQuest Mini, etc.).
2. Connect the interface to your device.
 - If using USB, connect to the USB port on your computer.
 - If using Bluetooth® wireless technology, click your interface type and then select your device.
3. Prepare for data collection:
 - Vernier Graphical Analysis®: Launch the app, if necessary, and click Sensor Data Collection.
 - LabQuest® App: Choose New from the File menu.

The software will identify the sensor and load a default data-collection setup. You are now ready to collect data.

⚠ WARNING: To avoid possible electric shock or personal injury, do not connect the red or black leads to household power. This product is designed to measure low-voltage sources such as classroom-scale wind turbines and small solar panels. It should never be connected to an electrical outlet.

Need Additional Information?

Visit the following link:

www.vernier.com/start-lq-sensor

Using the Product

Zeroing the Sensor

Both Current and Voltage should be zeroed in the data-collection software prior to collecting data. To do this, disconnect the load and source and connect the two

Load terminals with a wire. Zero the sensors in the software. **Note:** Any resistance value shown in the meter in the software is not meaningful when the current and voltage values are near zero.

Connecting Devices to the Energy Sensor

The Vernier Energy Sensor provides a pair of five-way binding posts, labeled Source, to connect the device that is generating power, such as a model wind turbine, solar panel, battery, power supply, or function generator. Another pair of binding post terminals, labeled Load, connect to the load. The load can be any electrical device that is meant to run on DC electricity at a voltage that matches the power source. Examples include the Vernier Variable Load (order code VES-VL), the Vernier Resistor Board (order code VES-RB), single component resistors, motors, or LEDs.

These terminals will accept standard banana plugs, alligator clips, or bare wires and provide a solid electrical connection.

Videos

View videos related to this product at www.vernier.com/ves-bta

Specifications

Source input potential range	± 30 V
Source input current range	± 1000 mA
Linearity	Potential Sensor: 0.01% Current Sensor: 0.01%
Resolution	Potential Sensor: 0.016 V Current Sensor: 0.52 mA
Input impedance	Potential Sensor: > 2 M Ω
Insertion resistance	Current Sensor: 0.1 Ω
Frequency response	Potential Sensor: -3dB @ 160 Hz Current Sensor: -3dB @ 160 Hz
Stored calibration (Potential Sensor)	slope: -12.526 V/V intercept: 31.315 V
Stored calibration (Current Sensor)	slope: -422.09 mA/V intercept: 1055.25 mA

How the Sensor Works


The Vernier Energy Sensor has two BTA connectors: a Voltage connector and a Current connector. If you connect only the Voltage connector or the Current connector, the Energy Sensor will function as a simple voltage or current sensor. When both connectors are connected to a Vernier interface, Graphical Analysis

and LabQuest app automatically identifies the sensor as the Energy Sensor and loads the appropriate data- collection setup for this sensor. Calculated columns and meters for power, resistance, and energy are created in the file.

If you are using other software listed in the Compatible Software section, you may have to manually set up the calculated columns and meters for power, resistance, and energy.

If you are using other software listed in the Compatible Software section, you may have to manually set up the calculated columns and meters for power, resistance, and energy.

The Vernier Energy Sensor measures the potential across the load, as well as the current through the load.



Circuit diagram of how the Vernier Energy Sensor is wired internally

Care and Maintenance

Do not wrap the cable tightly around the sensor for storage. Repeatedly doing so can irreparably damage the wires and is not covered under warranty.

Troubleshooting

For troubleshooting and FAQs, see www.vernier.com/til/3181

Repair Information

If you have watched the related product video(s), followed the troubleshooting steps, and are still having trouble with your Vernier Energy Sensor, contact Vernier Technical Support at support@vernier.com or call 888-837-6437. Support specialists will work with you to determine if the unit needs to be sent in for repair. At that time, a Return Merchandise Authorization (RMA) number will be issued and instructions will be communicated on how to return the unit for repair.

Accessories/Replacements

Item	Order Code
Vernier Variable Load	VES-VL
Vernier Resistor Board	VES-RB
KidWind Advanced Wind Experiment Kit	KW-AWX
KidWind Basic Wind Experiment Kit	KW-BWX
KidWind MINI Wind Turbine	KW-MWT
KidWind 2V/400mA Solar Panel	KW-SP2V

Warranty

Warranty information for this product can be found on the Support tab at www.vernier.com/ves-bta

General warranty information can be found at www.vernier.com/warranty



Vernier Science Education
13979 SW Millikan Way • Beaverton, OR 97005-2886
Toll Free (888) 837-6437 • (503) 277-2299 • Fax (503) 277-2440
info@vernier.com • www.vernier.com

Rev. 8/22/2024

Vernier Graphical Analysis, LabQuest, LabQuest Mini, and other marks shown are our trademarks or registered trademarks in the United States.

All other marks not owned by us that appear herein are the property of their respective owners, who may or may not be affiliated with, connected to, or sponsored by us.

