

NSTA 2025  
Philadelphia, PA

# Found a Box of Vernier Sensors in Your Classroom? Start Here!

## Experiments

### Introduction to Data Collection

- Go Direct® Temperature Probe

### Boyle's Law

- Go Direct Gas Pressure Sensor

## Workshop Presenter

Josh Ence  
[support@vernier.com](mailto:support@vernier.com)



Vernier Science Education • 888-VERNIER (888-837-6437)  
[info@vernier.com](mailto:info@vernier.com) • [www.vernier.com](http://www.vernier.com)





# Introduction to Data Collection

Data collection is a very important part of science. Meteorologists collect weather data over time to keep an historical record and to help make forecasts. Oceanographers collect data on the salinity (saltiness) of seawater to study changing trends in our Earth's oceans. While data have been collected by hand for thousands of years, the technology to collect data electronically was developed in the 1950s. Only since the 1980s has this technology been available and accessible to schools.

This experiment was designed to introduce you to two of the most common modes of data collection that will be used in this class. Part I will guide you through collecting and analyzing data over time. A temperature probe will be used to record the temperature of water for 60 seconds at a rate of one sample per second. In Part II, you will collect data using a mode called Events with Entry. This style of data collection allows you to collect one point of data, and then enter in a corresponding value. In this part, the data collected will be the temperature of your hand and the value you enter will be your assigned group member number.

## OBJECTIVES

- Become familiar with Graphical Analysis data-collection app.
- Use Graphical Analysis and a temperature probe to make measurements.
- Analyze a graph of the data.
- Use this graph to make conclusions about the experiment.
- Determine the response time of a temperature probe.

## MATERIALS


Chromebook, computer, **or** mobile device  
Graphical Analysis app  
Go Direct Temperature  
two 250 mL beakers  
cold tap water  
hot tap water  
ice  
paper towels




*Figure 1*

## PROCEDURE

### Part I Time Graph

1. Place about 100 mL of tap water into a 250 mL beaker. Add two or three ice cubes.
2. Launch Graphical Analysis. Connect the Go Direct Temperature Probe to your Chromebook, computer, or mobile device.
3. Set up the data-collection mode.
  - a. Click or tap Mode to open Data Collection Settings.
  - b. Change Rate to 0.5 samples/s and End Collection to 60 s. Click or tap Done.
4. Place the temperature probe into the cold water and stir briefly. Then position the probe in the cold-water beaker as shown in Figure 1. **Note:** Make sure the beaker will not tip over from the weight of the probe.
5. Place about 150 mL of hot water into a second 250 mL beaker.
6. When everything is ready, click or tap Collect to start data collection. Do not stir or move the water.
7. When exactly 10 seconds have gone by, quickly move the probe to the beaker containing hot water and allow data collection to continue. Do not stir the water or move the probe during the remainder of the data collection period.
8. Data collection will stop automatically after 60 seconds.
9. Remove the probe from the beaker and dry it with a paper towel.
10. Determine the elapsed time when the highest temperature was reached.
  - a. When data collection is complete, a graph of temperature vs. time is displayed. Click or tap the graph to examine the data. **Note:** You can also adjust the Examine line by dragging the line.
  - b. Find the highest temperature.
  - c. Record this temperature (round to the nearest 0.1°C) and the time when it was first reached in the data table.
11. Sketch or export an image of your graph according to your teacher's instructions.
12. You can confirm the time when the highest temperature was reached by viewing the data table.
  - a. Click or tap View, , and turn on the Data Table.
  - b. Find the time when the highest temperature was first reached. Did you get the same time both ways?

**Part II Events with Entry**

13. Click or tap File, , and choose New Experiment. Click or tap Sensor Data Collection.
14. Set up the data-collection mode.
  - a. Click or tap Mode to open Data Collection Settings. Change Mode to Event Based.
  - b. Enter **Member** as the Event Name and leave the Units field blank. Click or tap Done.
15. Assign numbers to the members of your group by age with the oldest being number one. Record the names in the data table for Part II. Add more lines if needed.
16. Click or tap Collect to start data collection.



*Figure 2*

17. Measure the hand temperature of the first group member.
  - a. Group member number one should pick up the probe and hold its tip in the palm of his or her hand as shown in Figure 2.
  - b. Watch the live temperature readout. When the temperature stops rising, click or tap Keep.
  - c. You will be prompted to enter a number. Enter **1** as the student's group member number, then click or tap Keep Point. The temperature and group member number have been saved.
18. Cool the probe down by placing it in the cold water from Part I. Monitor the temperature on the screen and remove it from the water when the temperature reaches 25°C.
19. Pass the probe to the next group member.
20. Repeat Steps 17–19 until every group member has his or her hand temperature recorded, entering the correct group member number for each person.
21. Click or tap Stop to stop data collection.
22. Determine each person's hand temperature by using one of the methods described in Steps 10 and 12. Record the values in the data table.
23. Sketch or export an image of your graph according to your teacher's instructions.



# Boyle's Law: Pressure-Volume Relationship in Gases

The primary objective of this experiment is to determine the relationship between the pressure and volume of a confined gas. The gas we use will be air, and it will be confined in a syringe connected to a gas pressure sensor (see Figure 1). When the volume of the syringe is changed by moving the piston, a change occurs in the pressure exerted by the confined gas. This pressure change will be monitored using a gas pressure sensor. It is assumed that temperature will be constant throughout the experiment. Pressure and volume data pairs will be collected during this experiment and then analyzed. From the data and graph, you should be able to determine what kind of mathematical relationship exists between the pressure and volume of the confined gas. Historically, this relationship was first established by Robert Boyle in 1662 and has since been known as Boyle's law.




Figure 1

## OBJECTIVES

- Use a gas pressure sensor and a gas syringe to measure the pressure of an air sample at several different volumes.
- Determine the relationship between pressure and volume of the gas.
- Describe the relationship between gas pressure and volume in a mathematical equation.
- Use the results to predict the pressure at other volumes.

## MATERIALS


Chromebook, computer, **or** mobile device  
Graphical Analysis app  
Go Direct Gas Pressure  
20 mL gas syringe

## PROCEDURE

1. Prepare the data-collection equipment and an air sample for data collection.
  - a. Launch Graphical Analysis. Connect the Gas Pressure Sensor to your Chromebook, computer, or mobile device.
  - b. With the 20 mL syringe disconnected from the Gas Pressure Sensor, move the piston of the syringe until the front edge of the inside black ring (indicated by the arrow in Figure 1) is positioned at the 10.0 mL mark.
  - c. Attach the 20 mL syringe to the valve of the Gas Pressure Sensor.
2. Set up the data-collection mode.
  - a. Click or tap Mode to open Data Collection Settings. Change Mode to Event Based.
  - b. Enter **Volume** as the Event Name and **mL** as the Units. Click or tap Done.
3. To obtain the best data possible, you will need to correct the volume readings from the syringe. Look at the syringe; its scale reports its own internal volume. However, that volume is not the total volume of trapped air in your system since there is a little bit of space inside the pressure sensor.

To account for the extra volume in the system, you will need to add 0.8 mL to your syringe readings. For example, with a 5.0 mL syringe volume, the total volume would be 5.8 mL. It is this total volume that you will need for the analysis.


4. You are now ready to collect pressure and volume data. It is easiest if one person takes care of the gas syringe and another enters volumes.
  - a. Click or tap Collect to start data collection.
  - b. Move the piston so the front edge of the inside black ring (see Figure 2) is positioned at the 5.0 mL line on the syringe. Hold the piston firmly in this position until the pressure value displayed on the screen stabilizes.
  - c. Click or tap Keep and enter **5.8**, the gas volume (in mL). Remember, you are adding 0.8 mL to the volume of the syringe for the total volume. Click or tap Keep Point to store this pressure-volume data pair.



*Figure 2*


- d. Continue this procedure using syringe volumes of 10.0, 12.5, 15.0, 17.5, and 20.0 mL.
  - e. Click or tap Stop to stop data collection.
5. When data collection is complete, a graph of pressure vs. volume will be displayed. To examine the data pairs on the displayed graph, tap any data point. As you tap each data point, the pressure and volume values are displayed to the right of the graph. Record the pressure and volume data values in your data table.



6. Based on the graph of pressure *vs.* volume, decide what kind of mathematical relationship exists between these two variables, direct or inverse. To see if you made the right choice:
  - a. Click or tap Graph Options, , and choose Apply Curve Fit.
  - b. Select Power as the curve fit and Dismiss the Curve Fit box. The curve fit statistics are displayed for the equation in the form

$$y = ax^b$$


where  $x$  is volume,  $y$  is pressure,  $a$  is a proportionality constant, and  $b$  is the exponent of  $x$  (volume) in this equation. **Note:** The relationship between pressure and volume can be determined from the value and sign of the exponent,  $b$ .

- c. If you have correctly determined the mathematical relationship, the regression line should very nearly fit the points on the graph (that is, pass through or near the plotted points).
  - d. Rescale the axes on your graph by clicking or tapping Graph Options, . Choose Edit Graph Options and set the x-axis to display 0 to 25 mL and the y-axis to display 0 to 300 kPa. Dismiss the Graph Options box.
  - e. (optional) Export, download, or print the graph with the curve fit displayed.
7. With the best-fit curve still displayed, proceed directly to the Processing the Data section.

## DATA AND CALCULATIONS

Volume (mL)	Pressure (kPa)	Constant, $k$ ( $P / V$ or $P \cdot V$ )




## PROCESSING THE DATA

1. With the best-fit curve still displayed, click or tap Graph Options, , and turn on Interpolate. Dismiss the box and click the graph to interpolate. Move along the regression line until the volume value is 5.0 mL. Note the corresponding pressure value. Now move to the point where the volume value is doubled (10.0 mL). What does your data show happens to the pressure when the volume is *doubled*? Show the pressure values in your answer.

## Boyle's Law: Pressure-Volume Relationship in Gases

- Using the same technique as in Question 1, what does your data show happens to the pressure if the volume is *halved* from 20.0 mL to 10.0 mL? Show the pressure values in your answer.
- Using the same technique as in Question 1, what does your data show happens to the pressure if the volume is *tripled* from 5.0 mL to 15.0 mL? Show the pressure values in your answer.
- From your answers to the first three questions *and* the shape of the curve in the plot of pressure *vs.* volume, do you think the relationship between the pressure and volume of a confined gas is direct or inverse? Explain your answer.
- Based on your data, what would you expect the pressure to be if the volume of the syringe was increased to 40.0 mL? Explain or show work to support your answer.
- Based on your data, what would you expect the pressure to be if the volume of the syringe was decreased to 2.5 mL? Explain or show work to support your answer.
- What experimental factors are assumed to be constant in this experiment?
- One way to determine if a relationship is inverse or direct is to find a proportionality constant,  $k$ , from the data. If this relationship is direct,  $k = P/V$ . If it is inverse,  $k = P \cdot V$ . Based on your answer to Question 4, choose one of these formulas and calculate  $k$  for the seven ordered pairs in your data table (divide or multiply the  $P$  and  $V$  values). Show the answers in the third column of the Data and Calculations table.
- How *constant* were the values for  $k$  you obtained in Question 8? Good data may show some minor variation, but the values for  $k$  should be relatively constant.
- Using  $P$ ,  $V$ , and  $k$ , write an equation representing Boyle's law. Write a verbal statement that correctly expresses Boyle's law.

## EXTENSION

- To confirm that an inverse relationship exists between pressure and volume, a graph of pressure *vs.* *reciprocal of volume* ( $1/\text{volume}$ ) may also be plotted. To do this, it is necessary to create a new column of data, reciprocal of volume, based on your original volume data:
  - Click or tap More Options, , in the Volume column header in the table. Choose Add Calculated Column.
  - Enter **1/volume** as the Name and **1/mL** as the Units.
  - Click or tap Insert Expression and choose A/X as the expression.
  - Enter **1** as Parameter A and select Volume as the Column.
  - Click or tap Apply.
- Plot a best-fit regression line on your graph of pressure *vs.*  $1/\text{volume}$ :
  - Click or tap Graph Options, , and choose Edit Graph Options.
  - Enter **0** as the value for both the Left value for the x-axis and the Bottom value for the y-axis.
  - Dismiss the Graph Options box. Your graph should now include the origin (0,0).
  - Click or tap Graph Options, , and choose Apply Curve Fit.

### ***Boyle's Law: Pressure-Volume Relationship in Gases***

- e. Select Linear as the curve fit and Dismiss the Curve Fit box. The linear-regression statistics are displayed in the form:

$$y = mx + b$$

where  $x$  is  $1/\text{volume}$ ,  $y$  is pressure,  $m$  is a proportionality constant, and  $b$  is the  $y$ -intercept.

- f. If the relationship between  $P$  and  $V$  is an inverse relationship, the graph of pressure vs.  $1/\text{volume}$  should be direct; that is, the curve should be linear and pass through (or near) the origin. Examine your graph to see if this is true for your data.