

Heart Rate Hand Grips

Order Code HR-GRIP



The Heart Rate Hand Grips require a Polar® Transmitter Module (HR-TRANS, not included). This combination is ideal for continuously monitoring heart rate before, during, and after exercise or while a person is stationary. Data are wirelessly transmitted to iPad®, LabQuest 2, or other mobile devices that are Bluetooth® Smart Ready. For information on how to use the Heart Rate Hand Grip with a Polar Transmitter Module please download the Go Wireless® Heart Rate user guide at www.vernier.com/gw-hr.

NOTE: Vernier products are designed for educational use. Our products are not designed nor are they recommended for any industrial, medical, or commercial process such as life support, patient diagnosis, control of a manufacturing process, or industrial testing of any kind.

Optional Accessories

Exercise Heart Rate Strap (order code: HR-STRAP)

A chest strap that is used with a Polar Transmitter Module (sold separately). This is a hands-free option for continuously monitoring heart rate. For more information, see www.vernier.com/hr-strap.

Polar Transmitter Module (order code: HR-TRANS)

A module that transmits heart rate data wirelessly to Bluetooth Smart Ready devices or devices equipped to receive 5 kHz RF signals. For more information, see www.vernier.com/hr-trans.

Heart Rate Receiver (order code: HR-REC)

A 5 kHz RF receiver attached to a BTA connector. This allows the Polar Transmitter Module to communicate with Vernier interfaces that do not support Bluetooth Smart Ready devices. For more information, see www.vernier.com/hr-rec.

Disposal Instruction


When disposing of this electronic product, do not treat it as household waste. Its disposal is subject to regulations that vary by country and region. This item should be given to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring that this product is disposed of correctly, you help prevent potential negative consequences on human health or on the environment. The recycling of materials will help to conserve natural resources. For more detailed information about recycling this product, contact your local city office or your disposal service.

The symbol, shown here, indicates that this product must not be disposed of in a standard waste container.



Warranty

Vernier warrants the Heart Rate Hand Grips to be free from defects in materials and workmanship for a period of five years from the date of shipment to the customer. This warranty does not cover damage to the product caused by abuse or improper use.



Measure. Analyze. Learn.™
Vernier Software & Technology

13979 S. W. Millikan Way • Beaverton, OR 97005-2886
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440
info@vernier.com • www.vernier.com

Rev. 11/19/2020

Go Wireless is our trademark in the United States. All other marks not owned by us that appear herein are the property of their respective owners, who may or may not be affiliated with, connected to, or sponsored by us.



Printed on recycled paper.