

Motion Encoder Receiver (MEC-BTD)

The Motion Encoder Receiver is a replacement part for the Motion Encoder System* (order code VDS-EC). The Receiver cannot be used on its own.

See www.vernier.com/vds-ec for a description and user manual for the Motion Encoder System.



WARNING: The Motion Encoder Receiver is not compatible with the Texas Instruments TI-Nspire Lab Cradle. Connecting the Receiver to the Lab Cradle will render the Cradle inoperative, requiring repair by Texas Instruments.



Measure. Analyze. Learn.™

Vernier Software & Technology

13979 S. W. Millikan Way • Beaverton, OR 97005-2886
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440
info@vernier.com • www.vernier.com

* Patent Pending

Rev. 10/8/2014